

How to prepare S-26 GOLD COMFORT®:



Always wash and dry hands before preparing formula.



Wash bottles, teats, caps, cups, knife and spoon.



After washing, sterilise items by boiling for 5 minutes. Keep bottles and cups capped until use.



Boil fresh drinking water for 5 minutes. Allow to cool.



Pour cooled boiled water into feeding bottle or cup.



Fill can's enclosed scoop with S-26 GOLD COMFORT powder and level off with the sterilised knife. Always add 1 scoop (8.7g) of powder for every 60mL of water.



Cap bottle and shake or stir formula in cup. Test temperature of milk on wrist before feeding. Formula left in the bottle or cup after a feed must be discarded.

- Prepare each bottle or cup individually.
- Consult feeding guide for quantities of water and formula.
- Ideally, formula should be prepared just prior to feeding. Otherwise, refrigerate prepared formula and use within 24 hours.